

Type MConnect



FOR THE COMMISSIONING AND DIAGNOSIS OF TFM AND TPM

Software for the configuration and diagnosis of monitoring systems TFM or TPM

- Display of the values being monitored
- Menu-based navigation
- Selection of values to be monitored and of alarms
- Configuration of additional functions
- Selection of dialogue language and of unit of measure for the volume flow rate

Application

- MConnect configuration software for use with monitoring devices TFM-1, TFM-2 and TPM
- PC software, available for download free of charge, to be used by customers
- Used for the configuration of monitoring devices for fume cupboards or rooms, monitoring of volume flow rate, face velocity, and differential pressure
- Functional testing and fault diagnosis
- Rapid documentation and creation of databases
- Upload/download and saving of data records

Special features

- Device configuration with graphical support
- Display of the operating values being monitored and of the actual values
- Monitoring of two independent values; the corresponding types of alarm can be set
- Special functions such as switch contacts for alarm suppression
- Creating backup files to build up a project database
- Selection of dialogue language, selection of unit of measure for volume flow rate

TROX Belgium N.V.



Paepsem Business Park
Boulevard Paepsem 18G
1070 Bruxelles
Tel: +0032 (0)2/522.07.80
Fax: +0032 (0)2/520.21.78

Online-Services

- > [TROX Academy](#)

- > [Votre contact](#)

- > [Formulaire de contact en ligne](#)

- > [Recherche Produit A-Z](#)

- > [Formulaire de demande du Tarif](#)

Service-Hotlines

TROX Belgium
+32 (0) 2/522.07.80

[Contact](#)